



Proenkephalin A 119-159 (penKid)

A unique biomarker for real-time assessment of kidney function

Table of Content

Abstract	5
Introduction	6
PenKid physiology in the kidney	7
Sphingotest® penKid® diagnostics	8
Published clinical indications	9
Kidney function	10
Subclinical AKI	11
Critically ill patients	12
Renal replacement therapy	13
Summary	14
Abbreviations	15
References	16
Disclaimer	19



Author:

Cristiane Carvalho-Wodarz, PhD, MBA
Medical Science Liaison Manager

cristiane.carvalho@sphingotec.com





Abstract

Acute kidney injury (AKI) affects many critically ill patients, making timely information on kidney function important for early intervention, such as renal replacement therapy (RRT) and nephroprotective strategies (1). The current AKI diagnostic criteria rely on the changes in serum creatinine and urine output, although they are neither sensitive nor specific for AKI (2).

An emerging body of evidence demonstrates that Proenkephalin A 119-159 (**penKid**), **a sensitive biomarker for kidney glomerular function**, overcomes these limitations (3). As a new biomarker for diagnosing AKI, penKid measured at admission to the emergency department (ED) (4, 5) or intensive care unit (ICU) (6, 7), enables timely anticipation of worsening renal function within 48 hours before current diagnostic criteria based on creatinine are met (8, 9).

More than 60 peer-reviewed publications and additional submitted clinical studies investigated the significance of penKid in more than 80,000 patients allowing a clear understanding of the value of using penKid in critical care medicine. Above all, penKid levels demonstrate a strong correlation with glomerular filtration rate (GFR) measured through iohexol/iothalamate clearance, providing real-time assessment of kidney function and offering a fast blood-based alternative for the in vivo measurement of true GFR (3, 10, 11).

In summary, penKid rapidly provides physicians at the ED and ICU with urgently needed information on kidney function, which is complementary to the current diagnostic toolbox. Thus, as studies show, penKid may aid in guiding nephrotoxic drug administration (12), contribute to clinical decisions around renal replacement therapy (RRT) by monitoring renal function under therapy (13), management of patients after major vascular surgery (14).

As CE₀₁₂₃ IVD marked quantitative immunoassay, **sphingotest® penKid® aids in the diagnosis of acute kidney injury in patients with sepsis or septic shock (15).**



Introduction

Impaired kidney function affects a broad number of critically ill patients and contributes to morbidity and mortality. In fact, one in three ICU patients develops AKI (16, 17), making it a global burden with estimated 13 million cases per year (16). AKI accounts for many complications in patients with sepsis or septic shock (1). It is furthermore concerning that around 30% of hospitalized AKI patients discharged from the hospital continue to experience with persisting renal dysfunction (18).

The underlying causes of AKI include, among others, sepsis, complex surgical interventions, cardiogenic shock, and acute heart failure (2). It is associated with a high mortality rate and extensive hospital costs. In the United States, additional expenses related to AKI complications are estimated to be approximately \$42,000 per case (19). Therefore, timely information on renal function is essential to initiate and adapt treatment options like RRT and nephroprotective strategies early. Although several additional biomarkers of kidney damage have been developed, they have not made it to clinical routine.

In clinical practice, serum creatinine (sCr) and urine output are the cornerstones of the current diagnostic approach for defining AKI (2). However, they give limited and delayed information on changes during kidney injury and have low sensitivity and specificity (2, 20).

Scientific experts agree on the significant need for new biomarkers that timely mirror kidney function, thereby improving the prediction and monitoring of AKI, major adverse kidney events (MAKE), and the necessity of RRT (20).

Translated into clinical practice these insights might provide benefit in the guidance to start or stop RRT, the guidance of nephrotoxic drug administration, hospitalization decisions after catheterization-laboratory procedures, the prediction of delayed graft function, and assessment of contrast-induced nephropathy. GFR is considered to be the best indicator of kidney function (21). However, the gold standard technique for determining the true GFR is based on a complex and resource-intensive procedure using inulin, iothalamate, or iohexol that are not feasible in daily practice. Today's routine standard relies on determining an estimated GFR (22).

This method is often based on sCr, which is known to be affected by non-renal factors (age, gender, muscle mass, and others) and therefore inaccurate. sCr is unable to detect mild kidney insufficiency because its levels only begin to rise above the normal value when approximately 50% of renal function is already impaired; this is known as the creatinine-blind area (23).


Scientific evidence showed that penKid predicts the future change in sCr up to two days in advance independently from inflammation or common comorbidities (e.g., hypertension, Diabetes mellitus) (3, 24).

PenKid is a stable prohormone fragment of the Enkephalin family, with a long in vivo half-life, is stable after collection and is not influenced by sex or age (10).

Studies have shown that penKid strongly correlates with the kidney function and measured GFR. As penKid does not bind to proteins in the plasma and it is freely filtered in the glomerulus, it is a promising biomarker for assessing kidney function in critically ill patients (3, 10, 24, 25).

PenKid may, therefore, provide physicians with urgently needed information on top of the standard of care, with a potential role to predict kidney recovery (6, 26).


In the consensus statement of the ADQI workgroup, penKid was recognized as a relevant functional biomarker able to detect subclinical AKI in critically ill patients (28).




PenKid physiology in the kidney

The enkephalins are endogenous opioids widely expressed in the human body with primary activity on delta opioid receptors (29). Interestingly, a high concentration of opioid receptors is found in the kidney (3, 30). Studies suggest a possible regulatory role for enkephalins in the kidneys, such as the induction of diuresis and natriuresis through receptor agonism (31) or inhibiting antidiuretic hormone activity (32). Several studies have shown a strong correlation of plasma penKid with the mGFR, which indicate that penKid might be a suitable and accurate maker to estimate the true GFR (25).

PenKid in the healthy state



PenKid in the disease state




Sphingotest[®] penKid[®] diagnostics



Proenkephalin A 119-159 is detected in the plasma with the diagnostic assay sphingotest[®] penKid[®]. The clinical performance of the CE₀₁₂₃ certified IVD sphingotest[®] penKid[®] was demonstrated in the patient cohort from the AdrenOSS-1 trial (i.e., adult patients with severe sepsis or septic shock) (35), using a cut-off at 89 pmol/L, resembling the upper normal limit of healthy subjects. Primary endpoint of this study was diagnosis of all-stage AKI within 48 hours after enrolment (n=529 patients).



Intended use: The sphingotest[®] penKid[®] is a non-automated immunoluminometric assay (ILMA) for the in vitro diagnostic quantitative measurement of Proenkephalin A 119-159 in human EDTA plasma. Proenkephalin A 119-159 measurement in conjunction with clinical assessments and other laboratory findings is used as an aid in the diagnosis of acute kidney injury in adult patients with sepsis or septic shock. The sphingotest[®] penKid[®] is intended for use by professional users in laboratory settings.

Sphingotest[®] penKid[®] performance in diagnosing AKI within 48 hours



Parameter	Observed performance
Diagnostic sensitivity	71.8% (95% CI 65.7%-77.2%)
Diagnostic specificity	83.1% (95% CI 78.4%-86.9%)
PPV	77.1% (95% CI 71.0%- 82.2%)
NPV	78.8% (95% CI 73.9%-83.0%)
Diagnostic odds ratio	12.5 (95% CI 8.2-18.9)

Figure 2. Sensitivity, specificity, PPV, NPV and diagnostic odds ratio for sphingotest[®] penKid[®] values for diagnosis of AKI within 48 hours. Clinical cut-off: 89 pmol/L.





Published clinical indications

Over 60 peer-reviewed publications have investigated penKid in more than 80,000 patients allowing a clear understanding of the characteristics of penKid levels and its value in critical care. PenKid has already compiled scientific evidence for many indications in critical care and beyond, proving value beyond standard of care markers.

Table 1: Main clinical indications where penKid showed additive value as a kidney function marker.

Investigated indications	PenKid value
Sepsis and septic shock	<p>ED: Earlier detection of AKI in sepsis patients upon ED admission is an unmet need for timely stratification of patients with AKI. PenKid showed great value by effectively detecting the presence and severity of AKI in septic patients with normal serum creatinine levels at admission (4, 5).</p> <p>ICU: Patients with sepsis or septic shock admitted to the ICU require efficient management strategies to predict kidney deterioration and recovery of renal function. PenKid predicted major adverse kidney events, transient AKI, worsening of renal function, and renal recovery in septic ICU patients (6, 7).</p>
Heart failure	<p>Acute heart failure (AHF): Patients with AHF, renal dysfunction, and inadequate decongestion during hospitalization are at risk of prolonged hospitalization, rehospitalization, and death. The rapid increase of penKid aids in identifying patients who do not tolerate the intensified diuretic therapy and predicts the worsening of renal function, in-hospital mortality, and mortality during follow-up (36-39).</p> <p>Chronic heart failure (CHF): There is an unmet clinical need for biomarkers to monitor kidney function in patients with CHF, mainly for patients with heart failure with preserved ejection fraction (HFpEF) (i.e., the left ventricle is not filling correctly with blood). High levels of penKid in HFpEF were associated with heart failure rehospitalization (40) and recurrence of acute myocardial infarction (41).</p> <p>Ambulatory patients: Patients with heart failure have a poor prognosis, but outcomes might be improved by early identification of organ failure risk. Increased penKid levels have shown significant predictive value for risk stratification of stable ambulatory patients (42).</p>
Major surgeries	<p>AKI occurs in up to 40% of patients undergoing cardiac surgery and is a common postoperative complication in complex aortic surgery. Due to its dynamic nature, penKid is a biomarker for the prediction of AKI before and after invasive interventions such as cardiac and aortic surgery (9, 14, 43).</p>
Organ transplantation	<p>PenKid is associated with kidney function in renal transplant recipients (RTR) and in healthy kidney donors and may be a marker for determining long-term renal prognosis in RTR. High penKid levels may indicate an additional benefit in the early identification of RTR at risk of graft failure (44). In liver transplantation during preoperative period, penKid was significantly higher in patients with severe AKI, and detected AKI after transplantation 48h earlier than sCr (8).</p>
RRT	<p>Critically ill patients are often at risk of developing AKI and are likely in need of organ support, such as RRT. However, the most effective timing to start RRT remains uncertain (45). PenKid levels have been shown to predict the need for RRT on day one or later during hospitalization in sepsis patients (7, 46) as well as in liver transplant patients (8). The ELAIN and RICH studies scientifically proved that penKid remains informational during RRT, indicating kidney function recovery and supporting the decision to safely stop RRT/CRRT (13, 47).</p>
Burned patients	<p>PenKid levels in patients on admission to the ICU were associated with the risk of developing AKI in patients with severe burns, where mortality rates range from 30-70%. The data suggest that penKid levels could provide additional value to the currently available markers of renal function and AKI to better quantify kidney function in critically ill patients (48).</p>
Pediatric ICU	<p>Pediatric AKI is often diagnosed too late for successful therapeutic interventions, leading to adverse outcomes and chronic kidney impairment in the long term (49). PenKid identifies AKI in critically ill neonates and children under one year of age, offering a potential solution to the shortcomings in pediatric critical care (50, 51).</p>



Kidney function

SCIENTIFIC EVIDENCE

PenKid has been evaluated as a biomarker of kidney function in several clinical conditions (6, 10, 52), and correlates with the gold standard measurement of GFR (10, 24). A penKid-based formula to estimate the GFR was developed and validated in a broad patient cohort (11). The assessment of penKid identifies patients with impaired kidney function, at high risk of AKI and delivers information on the GFR in healthy subjects and hospitalized patients (3, 6, 7). PenKid values reflect the measure and the true GFR and are independent of inflammation (7, 10, 11).

PenKid as a specific biomarker for glomerular filtration rate

PenKid is a reliable surrogate for the measured GFR in adults with stable and unstable renal conditions (10)

PenKid has a high correlation with the measured GFR in neonates and children (51)




Figure 5. PenKid inversely correlates with the measured GFR (iothalamate clearance): the higher the penKid levels, the lower the mGFR. N = 1191 patients. Modified from (10).




Figure 6. Correlation of penKid concentration as a function of measured GFR (mGFR, iothalamate clearance). N= 37 patients: $p = -0.88$, $p < 0.001$. Modified from (51).

PenKid eGFR formula correlates with mGFR (11)

PenKid has high specificity for renal function (7)




Figure 7. Correlation of mGFR and eGFR using the new penKid formula including penKid and sCr, developed from patients with corresponding assessment of GFR by iothexol or iothalamate clearance. Patient cohort: chronic kidney disease patients, healthy volunteers, sepsis and septic shock patients, kidney donors, and recipients. N=1354 patients. R-squared: 0.73 ($p < 0.001$). Modified from (11).




Figure 8. Association of penKid with renal SOFA score in severe septic and septic shock patients. PenKid increased only upon renal dysfunction (subscores 1 – 4), remaining within normal range even in the context of systemic inflammation, as in sepsis or septic shock (subscores 0). The line indicates the upper normal range. N = 956 patients. Modified from (7).



Subclinical AKI

SCIENTIFIC EVIDENCE

Subclinical AKI (sub-AKI) is defined as a condition where there is an increase of functional or damage biomarker, without signs of AKI according to the KDIGO criteria, i.e., sCr and urine output are still within normal range. The prevalence of sub-AKI on ICU admission is still not properly explored (between 12-18% patients) (28), however, there is evidence for poor prognosis in such patients. The current scientific evidence shows that penKid can detect sub-AKI, allowing early identification of patients at risk while not meeting current AKI criteria (8, 28, 48, 53).

PenKid identifies patients with sub-AKI at increased risk of death

PenKid identifies sub-clinical AKI in critically ill patients (28)




Figure 9. Survival of critically ill patients, based on presence of AKI, sub-AKI or no AKI. Sub-AKI defined based on penKid > 80 pmol/L. N=2087 patients. Modified from (28).

PenKid identifies sub-clinical AKI in burned patients (48)




Figure 10. Survival of burned patients, based on AKI, sub-AKI or no AKI. The identification of sub-AKI allows the identification of patients who do not fall into the clinical AKI definition but still show a GFR decrease and increased risk of death. Sub-AKI defined based on penKid > 80 pmol/L. N=113 patients. Modified from (48).




Critically ill patients


SCIENTIFIC EVIDENCE

The assessment of penKid identifies patients at high risk of AKI in EDs and ICUs and delivers information on kidney function in hospitalized patients (3-7, 54, 55).


PenKid can detect AKI and its severity in sepsis patients at the ED (4)




PenKid rises progressively with the severity of AKI in sepsis patients on the ICU (7)



High penKid levels on ED admission identify sepsis patients with poor outcome (54)



PenKid to monitor treatment success in septic patients (7)






Renal replacement therapy


SCIENTIFIC EVIDENCE

Renal replacement therapy (RRT) is the main therapy to treat AKI complications, such as electrolyte or acid-base imbalances, uremia, or fluid overload (13). Yet, the best time to start or stop RRT remains controversial (2), and the decision to stop RRT relies on the clinical assessment of patients, measurement of creatinine clearance and 24-h urine output – both not reliable parameters (56). PenKid has shown value in identifying patients requiring RRT and aiding in their successful RRT liberation (7, 13).

PenKid is higher in patients requiring RRT (7)




PenKid is higher in patients after liver transplantation receiving RRT (8)




PenKid is associated with the successful liberation* from RRT (13)

PenKid assessment in patients during the third day of RRT




Patients with increased penKid levels might benefit from early initiation of RRT



*Successful liberation: Patients who survived and did not receive any form of RRT for at least 7 days (relapse-free period) after CRRT discontinuation were classified as successfully liberated.

Summary

- PenKid is an innovative biomarker for the early diagnosis of AKI in critically ill patients.
- Publications show that penKid correlates with true GFR and enables early assessment of worsening and improving kidney function.
- PenKid has been measured in over 80,000 patients.
- Scientific evidence shows that penKid blood level is not influenced by inflammation.



Scheme displaying the rise of penKid blood levels (1) predicting worsening of kidney function up to 48 hours earlier than today's standard of care; the decrease of penKid levels (2) indicates the normalization of renal function.

Literature on penKid as a new functional kidney biomarker

Clinical settings already evaluated:

- Acute myocardial infarction (41)
- Burns (48)
- Cardiac arrest (57)
- Cardiac surgery (9, 14, 43, 58)
- Cardiogenic shock (59)
- Chronic kidney disease (60, 61)
- Contrast-induced nephrotoxicity (12)
- Heart failure (36-40, 42, 62-64)
- ICU (50, 65-67)
- Renal and liver transplantation (8, 44)
- RRT (7, 13, 47)
- Sepsis (4-7, 24, 28, 46, 52, 68, 69)

Other potential applications where penKid might benefit:

- Guidance of nephrotoxic drug administration.
- ICU discharge decisions.
- Management of patients attending catheterization intervention.
- Triaging.

Abbreviations

ADQI:	Acute Disease Quality Initiative
AHF:	Acute heart failure
AKI:	Acute kidney injury
CHF:	Chronic heart failure
ED:	Emergency department
EDTA:	Ethylenediaminetetraacetic acid
GFR:	Glomerular filtration rate
HFpEF:	Heart failure with preserved ejection fraction
ICU:	Intensive care unit
KDIGO:	Kidney disease: improving global outcomes
MAKE:	Major adverse kidney events
mGFR:	measured glomerular filtration rate
NPV:	Negative predictive value
PenKid:	Proenkephalin A 119-159
PPV:	Positive predictive value
RIFLE:	Risk, injury, failure, loss and end-stage kidney disease
ROC:	receiver operating characteristic
RRT:	Renal replacement therapy
RTR:	Renal transplants recipients
Sub-AKI:	Subclinical acute kidney injury
sCr:	Serum creatinine
SOFA:	Sequential Organ Failure Assessment

References

1. Bellomo R, Kellum JA, Ronco C, Wald R, Martensson J, Maiden M, et al. [Acute kidney injury in sepsis. Intensive Care Med. 2017;43\(6\):816-28.](#)
2. Ronco C, Bellomo R, Kellum JA. [Acute kidney injury. Lancet. 2019;394\(10212\):1949-64.](#)
3. Beunders R, Struck J, Wu AHB, Zarbock A, Di Somma S, Mehta RL, et al. [Proenkephalin \(PENK\) as a Novel Biomarker for Kidney Function. J Appl Lab Med. 2017;2\(3\):400-12.](#)
4. Marino R, Struck J, Hartmann O, Maisel AS, Rehfeldt M, Magrini L, et al. [Diagnostic and short-term prognostic utility of plasma pro-enkephalin \(pro-ENK\) for acute kidney injury in patients admitted with sepsis in the emergency department. J Nephrol. 2015;28\(6\):717-24.](#)
5. Rosenqvist M, Brnton K, Hartmann O, Bergmann A, Struck J, Melander O. [Proenkephalin a 119-159 \(penKid\) - a novel biomarker for acute kidney injury in sepsis: an observational study. BMC Emerg Med. 2019;19\(1\):75.](#)
6. Hollinger A, Wittebole X, François B, Pickkers P, Antonelli M, Gayat E, et al. [Proenkephalin A 119-159 \(Penkid\) Is an Early Biomarker of Septic Acute Kidney Injury: The Kidney in Sepsis and Septic Shock \(Kid-SSS\) Study. Kidney Int Rep. 2018;3\(6\):1424-33.](#)
7. Caironi P, Latini R, Struck J, Hartmann O, Bergmann A, Bellato V, et al. [Circulating Proenkephalin, Acute Kidney Injury, and Its Improvement in Patients with Severe Sepsis or Shock. Clin Chem. 2018;64\(9\):1361-9.](#)
8. Lima C, Gorab DL, Fernandes CR, Macedo E. [Role of proenkephalin in the diagnosis of severe and subclinical acute kidney injury during the perioperative period of liver transplantation. Pract Lab Med. 2022;31:e00278.](#)
9. Shah KS, Taub P, Patel M, Rehfeldt M, Struck J, Clopton P, et al. [Proenkephalin predicts acute kidney injury in cardiac surgery patients. Clin Nephrol. 2015;83\(1\):29-35.](#)
10. Donato LJ, Meeusen JW, Lieske JC, Bergmann D, Sparwaßer A, Jaffe AS. [Analytical performance of an immunoassay to measure proenkephalin. Clin Biochem. 2018;58:72-7.](#)
11. Beunders Rea. [ESICM LIVES 2021: Part 1. Intensive Care Medicine Experimental. 2021;9\(1\):51.](#)
12. Breidhardt T, Jaeger C, Christ A, Klima T, Mosimann T, Twerenbold R, et al. [Proenkephalin for the early detection of acute kidney injury in hospitalized patients with chronic kidney disease. Eur J Clin Invest. 2018;48\(10\):e12999.](#)
13. von Groote T, Albert F, Meersch M, Koch R, Porschen C, Hartmann O, et al. [Proenkephalin A 119-159 predicts early and successful liberation from renal replacement therapy in critically ill patients with acute kidney injury: a post hoc analysis of the ELAIN trial. Crit Care. 2022;26\(1\):333.](#)
14. Gombert A, Barbati M, Hartmann O, Schulte J, Simon T, Simon F. [Proenkephalin A 119-159 May Predict Post-operative Acute Kidney Injury and in Hospital Mortality Following Open or Endovascular Thoraco-abdominal Aortic Repair. Eur J Vasc Endovasc Surg. 2020;60\(3\):493-4.](#)
15. penKid(R) SR. Instructions for Use (IFU), IFU-PEK-001 (EUen). March 2022.
16. Ponce Dea. [Acute kidney injury: risk factors and management challenges in developing countries. Int J Nephrol Renovasc Dis. 2016;9:193-200.](#)
17. Hoste EA, Bagshaw SM, Bellomo R, Cely CM, Colman R, Cruz DN, et al. [Epidemiology of acute kidney injury in critically ill patients: the multinational AKI-EPI study. Intensive Care Med. 2015;41\(8\):1411-23.](#)
18. Kellum JA, Sileanu FE, Bihorac A, Hoste EA, Chawla LS. [Recovery after Acute Kidney Injury. Am J Respir Crit Care Med. 2017;195\(6\):784-91.](#)
19. Hall PS, Mitchell ED, Smith AF, Cairns DA, Messenger M, Hutchinson M, et al. [The future for diagnostic tests of acute kidney injury in critical care: evidence synthesis, care pathway analysis and research prioritisation. Health Technol Assess. 2018;22\(32\):1-274.](#)
20. Ostermann M, Zarbock A, Goldstein S, Kashani K, Macedo E, Murugan R, et al. [Recommendations on Acute Kidney Injury Biomarkers From the Acute Disease Quality Initiative Consensus Conference: A Consensus Statement. JAMA Netw Open. 2020;3\(10\):e2019209.](#)
21. Schaeffner E. [Determining the Glomerular Filtration Rate-An Overview. J Ren Nutr. 2017;27\(6\):375-80.](#)
22. Levey AS, Inker LA, Coresh J. [GFR estimation: from physiology to public health. Am J Kidney Dis. 2014;63\(5\):820-34.](#)
23. Shemesh O, Golbetz H, Kriss JP, Myers BD. [Limitations of creatinine as a filtration marker in glomerulopathic patients. Kidney Int. 1985;28\(5\):830-8.](#)
24. Beunders R, van Groenendaal R, Leijte GP, Kox M, Pickkers P. [Proenkephalin Compared to Conventional Methods to Assess Kidney Function in Critically Ill Sepsis Patients. Shock. 2020;54\(3\):308-14.](#)
25. Khorashadi M, Beunders R, Pickkers P, Legrand M. [Proenkephalin: A New Biomarker for Glomerular Filtration Rate and Acute Kidney Injury. Nephron. 2020;144\(12\):655-61.](#)
26. Patschan D, Erfurt S, Oess S, Lauxmann MA, Patschan S, Ritter O, et al. [Biomarker-based prediction of survival and recovery of kidney function in acute kidney injury. Kidney Blood Press Res. 2023.](#)
27. Zarbock A, Nadim MK, Pickkers P, Gomez H, Bell S, Joannidis M, et al. [Sepsis-associated acute kidney injury: consensus report of the 28th Acute Disease Quality Initiative workgroup. Nature Reviews Nephrology. 2023.](#)
28. Dépret F, Hollinger A, Cariou A, Deye N, Vieillard-Baron A, Fournier MC, et al. [Incidence and Outcome of Subclinical Acute Kidney Injury Using penKid in Critically Ill Patients. Am J Respir Crit Care Med. 2020;202\(6\):822-9.](#)
29. Grossman A, Clement-Jones V. [Opiate receptors: enkephalins and endorphins. Clin Endocrinol Metab. 1983;12\(1\):31-56.](#)

30. Denning GM, Ackermann LW, Barna TJ, Armstrong JG, Stoll LL, Weintraub NL, et al. [Proenkephalin expression and enkephalin release are widely observed in non-neuronal tissues.](#) *Peptides.* 2008;29(1):83-92.
31. Sezen SF, Kenigs VA, Kapusta DR. [Renal excretory responses produced by the delta opioid agonist, BW373U86, in conscious rats.](#) *J Pharmacol Exp Ther.* 1998;287(1):238-45.
32. Grossman A, Besser GM, Milles JJ, Baylis PH. [Inhibition of vasopressin release in man by an opiate peptide.](#) *Lancet.* 1980;2(8204):1108-10.
33. Ernst A, Köhrle J, Bergmann A. [Proenkephalin A 119-159, a stable proenkephalin A precursor fragment identified in human circulation.](#) *Peptides.* 2006;27(7):1835-40.
34. Fuchs MAA, Schrankl J, Wagner C, Daniel C, Kurtz A, Broeker KA. [Localization and characterization of proenkephalin-A as a potential biomarker for kidney disease in murine and human kidneys.](#) *Biomarkers.* 2023;28(1):76-86.
35. Mebazaa A, Geven C, Hollinger A, Wittebole X, Chousterman BG, Blet A, et al. [Circulating adrenomedullin estimates survival and reversibility of organ failure in sepsis: the prospective observational multinational Adrenomedullin and Outcome in Sepsis and Septic Shock-1 \(AdrenOSS-1\) study.](#) *Crit Care.* 2018;22(1):354.
36. Emmens JE, Ter Maaten JM, Damman K, van Veldhuisen DJ, de Boer RA, Struck J, et al. [Proenkephalin, an Opioid System Surrogate, as a Novel Comprehensive Renal Marker in Heart Failure.](#) *Circ Heart Fail.* 2019;12(5):e005544.
37. Molvin J, Jujic A, Navarin S, Melander O, Zoccoli G, Hartmann O, et al. [Bioactive adrenomedullin, proenkephalin A and clinical outcomes in an acute heart failure setting.](#) *Open Heart.* 2019;6(2):e001048.
38. Ng LL, Squire IB, Jones DJL, Cao TH, Chan DCS, Sandhu JK, et al. [Proenkephalin, Renal Dysfunction, and Prognosis in Patients With Acute Heart Failure: A GREAT Network Study.](#) *J Am Coll Cardiol.* 2017;69(1):56-69.
39. Matsue Y, Ter Maaten JM, Struck J, Metra M, O'Connor CM, Ponikowski P, et al. [Clinical Correlates and Prognostic Value of Proenkephalin in Acute and Chronic Heart Failure.](#) *J Card Fail.* 2017;23(3):231-9.
40. Kanagala P, Squire IB, Jones DJL, Cao TH, Chan DCS, McCann G, et al. [Proenkephalin and prognosis in heart failure with preserved ejection fraction: a GREAT network study.](#) *Clin Res Cardiol.* 2019;108(8):940-9.
41. Ng LL, Sandhu JK, Narayan H, Quinn PA, Squire IB, Davies JE, et al. [Proenkephalin and prognosis after acute myocardial infarction.](#) *J Am Coll Cardiol.* 2014;63(3):280-9.
42. Arbit B, Marston N, Shah K, Lee EL, Aramin H, Clopton P, et al. [Prognostic Usefulness of Proenkephalin in Stable Ambulatory Patients With Heart Failure.](#) *Am J Cardiol.* 2016;117(8):1310-4.
43. Mossanen JC, Pracht J, Jansen TU, Buendgens L, Stoppe C, Goetzenich A, et al. [Elevated Soluble Urokinase Plasminogen Activator Receptor and Proenkephalin Serum Levels Predict the Development of Acute Kidney Injury after Cardiac Surgery.](#) *Int J Mol Sci.* 2017;18(8).
44. Kieneker LM, Hartmann O, Struck J, Bergmann A, Gansevoort RT, Joosten MM, et al. [Plasma Proenkephalin and Poor Long-Term Outcome in Renal Transplant Recipients.](#) *Transplant Direct.* 2017;3(8):e190.
45. Boyer N, Horne K, Selby NM, Forni LG. [Renal medicine in the intensive care unit: a narrative review.](#) *Anaesthesia.* 2023;78(7):861-73.
46. Kim H, Hur M, Struck J, Bergmann A, Di Somma S. [Proenkephalin Predicts Organ Failure, Renal Replacement Therapy, and Mortality in Patients With Sepsis.](#) *Ann Lab Med.* 2020;40(6):466-73.
47. von Groote T, Albert F, Meersch M, Koch R, Gerss J, Arlt B, et al. [Evaluation of Proenkephalin A 119-159 for liberation from renal replacement therapy: an external, multicenter pilot study in critically ill patients with acute kidney injury.](#) *Crit Care.* 2023;27(1):276.
48. Dépret F, Polina A, Amzallag J, Fayolle-Pivot L, Coutrot M, Chaussard M, et al. [PenKid measurement at admission is associated with outcome in severely ill burn patients.](#) *Burns.* 2020;46(6):1302-9.
49. Lebel A, Teoh CW, Zappitelli M. [Long-term complications of acute kidney injury in children.](#) *Curr Opin Pediatr.* 2020;32(3):367-75.
50. Hartman SJF, Zwiers AJM, van de Water NEC, van Rosmalen J, Struck J, Schulte J, et al. [Proenkephalin as a new biomarker for pediatric acute kidney injury - reference values and performance in children under one year of age.](#) *Clin Chem Lab Med.* 2020;58(11):1911-9.
51. Smeets NJL, Hartmann O, Schulte J, Schreuder MF, de Wildt SN. [Proenkephalin A as a marker for glomerular filtration rate in critically ill children: validation against gold standard iohexol GFR measurements.](#) *Clin Chem Lab Med.* 2023;61(1):104-11.
52. Kim H, Hur M, Lee S, Marino R, Magrini L, Cardelli P, et al. [Proenkephalin, Neutrophil Gelatinase-Associated Lipocalin, and Estimated Glomerular Filtration Rates in Patients With Sepsis.](#) *Ann Lab Med.* 2017;37(5):388-97.
53. Zou C, Wang C, Lu L. [Advances in the study of subclinical AKI biomarkers.](#) *Front Physiol.* 2022;13:960059.
54. Papisidero ID, Valli G, Marin D, Del Sasso A, De Magistris A, Cennamo E, et al. [Utility of Measuring Circulating Bio-Adrenomedullin and Proenkephalin for 30-Day Mortality Risk Prediction in Patients with COVID-19 and Non-COVID-19 Interstitial Pneumonia in the Emergency Department.](#) *Medicina (Kaunas).* 2022;58(12).

-
55. Casalbani S, Valli G, Terlizzi F, Mastracchi M, Fidelio G, De Marco F, et al. [30 Days Mortality Prognostic Value of POCT Bio-Adrenomedullin and Proenkephalin in Patients with Sepsis in the Emergency Department. *Medicina \(Kaunas\)*. 2022;58\(12\).](#)
 56. Schiffli H, Lang SM. [Current Approach to Successful Liberation from Renal Replacement Therapy in Critically Ill Patients with Severe Acute Kidney Injury: The Quest for Biomarkers Continues. *Molecular Diagnosis & Therapy*. 2021;25\(1\):1-8.](#)
 57. Thorgeirsdóttir B, Levin H, Spångfors M, Annborn M, Cronberg T, Nielsen N, et al. [Plasma proenkephalin A 119-159 and dipeptidyl peptidase 3 on admission after cardiac arrest help predict long-term neurological outcome. *Resuscitation*. 2021;163:108-15.](#)
 58. Hill A, Bergmann D, Schulte J, Zayat R, Marx G, Simon TP, et al. [Proenkephalin A and bioactive adrenomedullin are useful for risk prognostication in cardiac surgery. *Front Cardiovasc Med*. 2022;9:1017867.](#)
 59. Jäntti T, Tarvasmäki T, Harjola VP, Pulkki K, Turkia H, Sabell T, et al. [Predictive value of plasma proenkephalin and neutrophil gelatinase-associated lipocalin in acute kidney injury and mortality in cardiogenic shock. *Ann Intensive Care*. 2021;11\(1\):25.](#)
 60. Kieneker LM, Hartmann O, Bergmann A, de Boer RA, Gansevoort RT, Joosten MM, et al. [Proenkephalin and risk of developing chronic kidney disease: the Prevention of Renal and Vascular End-stage Disease study. *Biomarkers*. 2018;23\(5\):474-82.](#)
 61. Schulz CA, Christensson A, Ericson U, Almgren P, Hindy G, Nilsson PM, et al. [High Level of Fasting Plasma Proenkephalin-A Predicts Deterioration of Kidney Function and Incidence of CKD. *J Am Soc Nephrol*. 2017;28\(1\):291-303.](#)
 62. Arfsten H, Goliash G, Bartko PE, Prausmüller S, Spinka G, Cho A, et al. [Neprilysin inhibition does not alter dynamic of proenkephalin-A 119-159 and pro-substance P in heart failure. *ESC Heart Fail*. 2021;8\(3\):2016-24.](#)
 63. Siong Chan DC, Cao TH, Ng LL. [Proenkephalin in Heart Failure. *Heart Fail Clin*. 2018;14\(1\):1-11.](#)
 64. Wu AHB, Anand I. [The biological variation of plasma proenkephalin: data from a stable heart failure cohort. *Clin Chem Lab Med*. 2019;57\(6\):e105-e7.](#)
 65. Frigyesi A, Boström L, Lengquist M, Johnsson P, Lundberg OHM, Spångfors M, et al. [Plasma proenkephalin A 119-159 on intensive care unit admission is a predictor of organ failure and 30-day mortality. *Intensive care medicine experimental \[Internet\]*. 2021 2021/07//; 9\(1\):\[36 p.\]](#)
 66. Legrand M, Hollinger A, Vieillard-Baron A, Depret F, Cariou A, Deye N, et al. [One-Year Prognosis of Kidney Injury at Discharge From the ICU: A Multicenter Observational Study. *Crit Care Med*. 2019;47\(12\):e953-e61.](#)
 67. Gayat E, Touchard C, Hollinger A, Vieillard-Baron A, Mebazaa A, Legrand M, et al. [Back-to-back comparison of penKID with NephroCheck® to predict acute kidney injury at admission in intensive care unit: a brief report. *Critical Care*. 2018;22\(1\):24.](#)
 68. Liu R, Zheng X, Wang H, Wang S, Yu K, Wang C. [The value of plasma pro-enkephalin and adrenomedullin for the prediction of sepsis-associated acute kidney injury in critically ill patients. *Crit Care*. 2020;24\(1\):162.](#)
 69. Moledina DG. [Penkid: A Novel Biomarker of Reduced GFR in Sepsis. *Kidney Int Rep*. 2019;4\(1\):17-9.](#)

Disclaimer

This communication is for educational use only. It is not intended to guide or imply medical decision making.

sphingotest® penKid® is offered for in vitro diagnostics. “PenKid” represents the analyte Proenkephalin A 119-159.

Sphingotest® penKid® is CE₀₁₂₃ IVD marked and hence certified and approved for use in the European Economic Area (EEA) only. Human diagnostic use of such products may be subject to local regulations.

The information contained in this communication does not constitute nor imply an offer to sell or transfer any product based on the SphingoTec MTP (sphingotest®) as in vitro diagnostic (IVD) product in the United States of America or Canada. No product based on the SphingoTec MTP is currently available for sale as an IVD product in the United States of America or Canada. The analytical and clinical performance characteristics in compliance with U.S. FDA medical device regulations of any product that may be sold at some future point in time in the U.S., have not yet been established.

For more information, please visit www.sphingotec.com

SphingoTec GmbH
Neuendorfstr. 15A
16761 Hennigsdorf, Germany
Phone: +49 3302 20565 0
Email: info@sphingotec.com
www.sphingotec.com