



Proenkephalin A 119-159 (penKid)

A unique biomarker for real-time assessment of kidney function

Table of Content

Abstract	5
Introduction	6
PenKid physiology in the kidney	7
Sphingotest® penKid®	8
Kidney function	9
Moving beyond creatinine	10
Critically ill patients	11
Kidney replacement therapy	12
Kidney transplantation	13
Summary	14
Abbreviations	15
References	16
Disclaimer	19





Abstract

Acute kidney injury (AKI) affects many critically ill patients, making timely information on kidney function important for early intervention, such as kidney replacement therapy (KRT) and nephroprotective strategies (1). The current AKI diagnostic criteria rely on the changes in serum creatinine and urine output, although they are neither sensitive nor specific for AKI (2).

An emerging body of evidence demonstrates that Proenkephalin A 119-159 (penKid), a dynamic biomarker for kidney glomerular function, overcomes these limitations (3). As a biomarker for diagnosing AKI, penKid measured at admission to the emergency department (ED) (4, 5) or intensive care unit (ICU) (6-9), enables timely anticipation of worsening renal function within 48 hours before current diagnostic criteria based on creatinine are met (10, 11).

More than 70 peer-reviewed publications report on the significance of penKid in more than 80,000 patients, allowing a clear understanding of the value of using penKid in critical care medicine. Above all, penKid levels demonstrate a strong correlation with glomerular filtration rate (GFR) measured through iohexol/iothalamate clearance, providing real-time assessment of kidney function and offering a fast blood-based alternative for the in vivo measurement of true GFR (3, 12, 13).

In summary, penKid provides physicians at the ED and ICU with critical information on kidney function, complementing the current diagnostic tools. As studies have shown, penKid contributes to clinical decision-making for patients under KRT by enabling the monitoring of renal function under therapy (8, 14, 15), and supports the management of patients following kidney transplantation (16, 17).



Introduction

Impaired kidney function affects a broad number of critically ill patients and contributes to morbidity and mortality. In fact, one in three ICU patients develops AKI (18, 19), making it a global burden with estimated 13 million cases per year (18). AKI accounts for many complications in patients with sepsis or septic shock (1). It is furthermore concerning that approximately 30% of hospitalized patients with AKI discharged from the hospital continue to experience persisting renal dysfunction (20).

The underlying causes of AKI include, among others, sepsis, complex surgical interventions, cardiogenic shock, and acute heart failure (2). It is associated with a high mortality rate and extensive hospital costs. In the United States, additional expenses related to AKI complications are estimated to be approximately \$42,000 per case (21). Therefore, timely information on renal function is essential to initiate and adapt treatment options like KRT and nephroprotective strategies early. Although several additional biomarkers of kidney damage have been developed, they have not made it to clinical routine.

In clinical practice, serum creatinine (sCr) and urine output are the cornerstones of the current diagnostic approach for defining AKI (2). However, they give limited and delayed information on changes during kidney injury and have low sensitivity and specificity (2, 22).

Scientific experts agree on the significant need for new biomarkers that timely mirror kidney function, thereby improving the prediction and monitoring of AKI, major adverse kidney events (MAKE), graft failure and the necessity of KRT (22, 23). When translated into clinical practice, these insights may inform decision-making related to the discontinuation of KRT, the guidance of nephrotoxic drug administration, hospitalization decisions after catheterization-laboratory procedures, the prediction of delayed graft function in kidney transplant recipients, and assessment of contrast-induced nephropathy. GFR is considered to be the best indicator of kidney function (24).

However, the gold standard technique for determining the true GFR is based on a complex and resource-intensive procedure using inulin, iothalamate, or iohexol that are not feasible in daily practice. Today's routine standard relies on determining an estimated GFR (25).

This method is often based on sCr, which is known to be affected by non-renal factors (age, gender, muscle mass, and others) and therefore inaccurate. SCr is unable to detect mild kidney insufficiency because its levels only begin to rise above the normal value when approximately 50% of renal function is already impaired; this is known as the creatinine-blind area (26).


Scientific evidence showed that penKid predicts the future change in sCr up to two days in advance independently from inflammation or common comorbidities (e.g., hypertension, Diabetes mellitus) (3, 27).

PenKid is a stable prohormone fragment of the enkephalin family, with a long in vivo half-life. It is stable after collection and is not influenced by sex or age in adults (12).

Studies have shown that penKid strongly correlates with kidney function and measured GFR. PenKid does not bind to plasma proteins and is freely filtered in the glomerulus, making it a promising biomarker for assessing kidney function in critically ill patients (3, 12, 27, 28).

PenKid may, therefore, provide physicians with urgently needed information on top of the standard of care, with a potential role to predict kidney recovery (6, 29).


In the consensus statement of the ADQI workgroup, penKid was recognized as a relevant functional biomarker able to detect sub-clinical AKI in critically ill patients (30).



PenKid physiology in the kidney


Enkephalins are endogenous opioids widely expressed in the human body with primary activity on delta opioid receptors (31). Interestingly, a high concentration of opioid receptors is found in the kidney (3, 32). Studies suggest a possible regulatory role for enkephalins in the kidneys, such as the induction of diuresis and natriuresis through receptor agonism (33) or inhibition of antidiuretic hormone activity (34). Several studies have shown a strong correlation of plasma penKid with mGFR, indicating that penKid might be a suitable and accurate surrogate maker to estimate the true GFR (28).

PenKid in the healthy state



Kidney function is modulated by the hormone enkephalin which remains difficult to detect due to its short half-life (12). PenKid overcomes this limitation, as it is the stable fragment of the Proenkephalin A precursor protein, therefore being a suitable surrogate measure of enkephalins (35).

PenKid in the disease state



Impaired renal function triggers the expression of Proenkephalin A, which is causative for the increase in plasma penKid levels. Thus, high penKid levels indicate impaired kidney function (28, 35).

Sphingotest[®] penKid[®]



Research use only product developed by SphingoTec

The research use only sphingotest[®] penKid[®] is a microtiter plate assay designed for the measurement of Proenkephalin A 119-159 in human plasma. The 97.5th percentile in healthy individuals is 89 pmol/L (Figure 1).




Figure 1. PenKid reference range of 120 healthy adults based on European population; age 31-82 (mean 51) years, 46% female. Median 29.9 pmol/L; 97.5th percentile: 89 pmol/L (95% CI 85-118) (36).

Microtiter plate (MTP)-based immunoassay sphingotest[®] penKid[®] (12)




Figure 2. The MTP is an appropriate platform for manual high-throughput analysis for research use. Results are obtained in 1 hour. Detection range: 30 - 1,172 pmol/L. Sample volume: 50 µL EDTA plasma (12).

Diagnostic solutions of our licensing partners

SphingoTec collaborates with leading diagnostic companies to bring advanced biomarker assays to laboratories and clinics around the world. Through strategic licensing agreements, our partners make SphingoTec's biomarkers available on their proprietary diagnostic platforms. These assays combine SphingoTec's scientific innovation with the manufacturing and distribution capabilities of our partners, offering healthcare professionals reliable, widely accessible solutions for critical patient needs.



Figure 3. AFIAS sphingotest[®] penKid[®] assay on the AFIAS diagnostic platform. This CE IVDR-certified fluorescence immunoassay quantitatively measures Proenkephalin A 119-159 (penKid). Commercialized by Bodi-tech Med Inc., a global leader in point-of-care diagnostics.



Figure 4. IB10 sphingotest[®] penKid[®] assay on the Nexus IB10 diagnostic platform. This CE IVD-released immunoassay quantitatively measures Proenkephalin A 119-159 (penKid). Commercialized by Nexus Dx, a provider of near-patient testing systems and advanced diagnostic solutions.



Kidney function



SCIENTIFIC EVIDENCE

PenKid has been evaluated as a biomarker of kidney function in several clinical conditions (6, 12, 37), and correlates with the gold standard measurement of GFR (12, 27). A penKid-based formula to estimate the GFR was developed and validated in a broad patient cohort (13). The assessment of penKid identifies patients with impaired kidney function, at high risk of AKI and delivers information on the GFR in healthy subjects and hospitalized patients (3, 6, 7). PenKid values reflect the measured GFR and are independent of inflammation (7, 12, 13).


PenKid as a surrogate biomarker for glomerular filtration rate

PenKid is a reliable surrogate for the measured GFR in adults with stable and unstable renal conditions (12).


PenKid has a high correlation with the measured GFR in neonates and children (38).



PenKid eGFR formula correlates with mGFR (13)



PenKid has high specificity for renal function (7)





Moving beyond creatinine

SCIENTIFIC EVIDENCE

The current diagnostic criteria for AKI routinely employ serum creatinine and urine output (39). However, both criteria are known to have several limitations. An increase in serum creatinine occurs only after more than 50% of the GFR is lost (40), urging the need for novel biomarkers for reliable and timely detection of kidney dysfunction. Subclinical AKI (sub-AKI) is characterized by an increase in functional or damage biomarkers in the absence of clinical criteria of AKI. The prevalence of sub-AKI on ICU admission is still not properly explored (between 12-18% patients) (30), however existing evidence suggests an association with poor prognosis in affected patients. The current scientific evidence shows that penKid can detect sub-AKI, allowing early identification of patients at risk while not meeting current AKI criteria (10, 30, 40, 41).

PenKid identifies patients with sub-AKI and increased risk of death

PenKid identifies sub-AKI in critically ill patients (30)




Figure 9. Survival of critically ill patients, based on presence of AKI, sub-AKI or no AKI. Sub-AKI defined based on penKid > 80 pmol/L. N=2087 patients. Modified from (30).

PenKid identifies sub-AKI in burned patients (41)




Figure 10. Survival of burned patients, based on AKI, sub-AKI or no AKI. The identification of sub-AKI allows the identification of patients who do not fall into the clinical AKI definition but still show a GFR decrease and increased risk of death. Sub-AKI defined based on penKid > 80 pmol/L. N=113 patients. Modified from (41).




Critically ill patients


SCIENTIFIC EVIDENCE

The assessment of penKid identifies patients at high risk of AKI and delivers information on kidney function in hospitalized patients (3-7, 23, 42, 43).


PenKid can detect AKI and its severity in sepsis patients at the ICU (8)




PenKid rises progressively with the severity of AKI in sepsis patients on the ICU (7)



High penKid levels on ED admission identify sepsis patients with poor outcome (43)



PenKid to monitor treatment success in septic patients (7)





Kidney replacement therapy

SCIENTIFIC EVIDENCE

Kidney replacement therapy (KRT) is the main therapy to correct severe complications of AKI complications, such as electrolyte or acid-base imbalances, uremia, or fluid overload (14). The monitoring of kidney function under KRT is crucial to guide successful discontinuation of therapy. To date, the decision to discontinue KRT relies on the clinical assessment of patients, measurement of creatinine clearance and urine output – parameters which are not reliable for profound decision-making (44). PenKid is freely filtered by the glomerulus (45) and effectively cleared by extracorporeal therapies including convection (CVVH), diffusion (CVVHD) and adsorption (HA) (46). Despite this, plasma concentrations of patients receiving therapy are maintained, making it a valuable tool for monitoring patients under KRT and for guiding KRT discontinuation, thereby supporting successful liberation from KRT (7, 14).

PenKid is higher in patients requiring KRT (7)




Figure 15. PenKid is higher in severe sepsis or septic shock patients requiring immediate KRT on day 1 or any time later. ALBIOS study. N=656 patients, $p<0.0001$. Modified from (7).

PenKid is a reliable kidney marker under KRT (8)




Figure 16. PenKid and serum creatinine trajectories in patients receiving KRT. Serum creatinine declines during KRT, whereas penKid continues to increase in these patients. N=137 patients. Modified from (8)

PenKid is associated with the successful liberation* from KRT (15)




Figure 17. Receiver-operating characteristics analysis for penKid at CKRT discontinuation to determine successful CKRT discontinuation. N=20 patients (15).




Figure 18. PenKid level at CKRT discontinuation discriminating between patients successfully (dark blue) and unsuccessfully (light blue) liberated from CKRT. N=20 patients, $p<0.05$. Modified from (15).

*Successful liberation: Patients who survived and did not receive any form of KRT for at least 7 days (relapse-free period) after CKRT discontinuation were classified as successfully liberated.

Kidney transplantation

SCIENTIFIC EVIDENCE

Kidney transplantation remains the gold standard treatment for end-stage renal disease. However, early post-transplant management is critical, particularly when it comes to identifying patients at risk for delayed graft function (DGF), a common complication after kidney transplantation (47). PenKid provides a significant clinical advantage by identifying patients at risk for DGF well before traditional markers like serum creatinine show any changes (16). Moreover, penKid distinguishes between slow and delayed graft function with remarkable accuracy, up to eight days earlier than conventional assessments (16). This early and precise insight supports timely and personalized treatment decisions, giving clinicians a powerful tool to improve transplant outcomes.

PenKid identifies slow- and delayed graft function post transplant




Figure 19. PenKid trajectories in kidney transplant recipients stratified by recovery of graft function. IGF - immediate graft function, SGF - slow graft function, DGF - delayed graft function. N=159 patients. * $P \leq 0.05$, ** $P \leq 0.01$, *** $P \leq 0.001$. Modified from (16).

PenKid indicates severity of delayed graft function post transplant




Figure 20. PenKid trajectories in kidney transplant recipients stratified by delayed graft function severity. DGF – delayed graft function. N=159 patients. * $P \leq 0.05$, ** $P \leq 0.01$, *** $P \leq 0.001$. Modified from (16).

PenKid predicts delayed graft function




Figure 21. Receiver-operating characteristics analysis to predict delayed graft function. Calculated is the relative change in penKid and serum creatinine from pre-transplant to first post-transplant day. N=159 patients receiving a kidney transplantation. Modified from (16).

High penKid levels indicate a higher risk of graft failure




Figure 22. Kaplan-Meier curve for graft failure per tertiles of penKid. Black line: penKid >134 pmol/L, grey line: penKid <38 pmol/L, dotted line: penKid 38-134 pmol/L. N=664 kidney transplant recipients. Modified from (17).

****Mild DGF:** Patients requiring kidney replacement therapy (KRT) within the first 24 h; **Moderate DGF:** Patients requiring KRT up to Day 7 post-transplant; **Severe DGF:** Patients requiring KRT beyond Day 7 post-transplant.



Summary

Cutting edge research highlights on penKid

- PenKid is an innovative biomarker for the early diagnosis of AKI in critically ill patients.
- PenKid correlates with true GFR.
 - PenKid is reported to assess renal function under KRT.
 - PenKid indicates delayed graft function on day 1 post-transplant.
- PenKid is not influenced by inflammation in critically ill patients including sepsis and septic shock.
- Scientific evidence on penKid relies on measurement in over 80,000 patients.

Table 1: Pooled performance criteria of penKid for early detection of AKI. Results from a systematic review with meta-analysis incorporating 11 observational studies with 3,969 patients and 23.4% AKI incidence (48).

Pooled sensitivity	0.69 (95% CI 0.62-0.75)
Pooled specificity	0.76 (95% CI 0.68-0.82)
Positive likelihood ratio	2.83 (95% CI 2.06-3.88)
Negative likelihood ratio	0.41 (95% CI 0.33-0.52)
Area under the curve	0.77 (95% CI 0.73-0.81)



Scheme displaying the rise of penKid blood levels [1] predicting worsening of kidney function up to 48 hours earlier than today's standard of care; the decrease of penKid levels [2] indicates the normalization of renal function.

Abbreviations

ADQI:	Acute Disease Quality Initiative
AKI:	Acute kidney injury
CKRT:	Continuous kidney replacement therapy
DGF:	Delayed graft function
ED:	Emergency department
GFR:	Glomerular filtration rate
ICU:	Intensive care unit
IGF:	Immediate graft function
KDIGO:	Kidney disease: improving global outcomes
KRT:	Kidney replacement therapy
MAKE:	Major adverse kidney events
mGFR:	measured glomerular filtration rate
NPV:	Negative predictive value
penKid:	Proenkephalin A 119-159
PPV:	Positive predictive value
ROC:	receiver operating characteristic
SA-AKI:	Sepsis associated- acute kidney injury
sCr:	Serum creatinine
SGF:	Slow graft function
SOFA:	Sepsis-relate organ failure assessment
WRF:	Worsening of renal function

References

1. Bellomo R, Kellum JA, Ronco C, Wald R, Martensson J, Maiden M, et al. Acute kidney injury in sepsis. *Intensive care medicine*. 2017;43(6):816-28.
2. Ronco C, Bellomo R, Kellum JA. Acute kidney injury. *The Lancet*. 2019;394(10212):1949-64.
3. Beunders R, Struck J, Wu AH, Zarbock A, Di Somma S, Mehta RL, et al. Proenkephalin (PENK) as a novel biomarker for kidney function. *J Appl Lab Med*. 2017;DOI: 10.1373/jalm.2017.023598 400-12.
4. Marino R, Struck J, Hartmann O, Maisel AS, Rehfeldt M, Magrini L, et al. Diagnostic and short-term prognostic utility of plasma pro-enkephalin (pro-ENK) for acute kidney injury in patients admitted with sepsis in the emergency department. *Journal of nephrology*. 2015;28(6):717-24.
5. Rosenqvist M, Bronton K, Hartmann O, Bergmann A, Struck J, Melander O. Proenkephalin A 119-159 (penKid) - a novel biomarker for acute kidney injury in sepsis: an observational study. *BMC Emerg Med*. 2019;19(1):75.
6. Hollinger A, Wittebole X, François B, Pickkers P, Antonelli M, Gayat E, et al. Proenkephalin A 119-159 (Penkid) Is an Early Biomarker of Septic Acute Kidney Injury: The Kidney in Sepsis and Septic Shock (Kid-SSS) Study. *Kidney Int Rep*. 2018;3(6):1424-33.
7. Caironi P, Latini R, Struck J, Hartmann O, Bergmann A, Bellato V, et al. Circulating Proenkephalin, Acute Kidney Injury, and Its Improvement in Patients with Severe Sepsis or Shock. *Clinical chemistry*. 2018;64(9):1361-9.
8. Martin L, Martin C, Peine A, Imohl M, Kersten A, Kramann R, et al. Implementation and One-Year Evaluation of Proenkephalin A in Critical Care. *International journal of molecular sciences*. 2025;26(6).
9. Schulte J, Dépret F, Hartmann O, Pickkers P, Laterre P-F, Uhle F, et al. Clinical performance of proenkephalin A 119-159 for the early diagnosis of acute kidney injury in patients with sepsis or septic shock. *medRxiv preprint*. 2024.
10. Lima C, Gorab DL, Fernandes CR, Macedo E. Role of proenkephalin in the diagnosis of severe and subclinical acute kidney injury during the perioperative period of liver transplantation. *Pract Lab Med*. 2022;31:e00278.
11. Shah KS, Taub P, Patel M, Rehfeldt M, Struck J, Clopton P, et al. Proenkephalin predicts acute kidney injury in cardiac surgery patients. *Clinical nephrology*. 2015;83(1):29-35.
12. Donato LJ, Meeusen JW, Lieske JC, Bergmann D, Sparwaßer A, Jaffe AS. Analytical performance of an immunoassay to measure proenkephalin. *Clin Biochem*. 2018;58:72-7.
13. Beunders R, Donato LJ, van Groenendael R, Arlt B, Carvalho-Wodarz C, Schulte J, et al. Assessing GFR With Proenkephalin. *Kidney Int Rep*. 2023;8(11):2345-55.
14. von Groote T, Albert F, Meersch M, Koch R, Porschen C, Hartmann O, et al. Proenkephalin A 119-159 predicts early and successful liberation from renal replacement therapy in critically ill patients with acute kidney injury: a post hoc analysis of the ELAIN trial. *Critical care*. 2022;26(1):333.
15. Tichy J, Hausmann A, Lanzerstorfer J, Ryz S, Wagner L, Lassnigg A, et al. Prediction of Successful Liberation from Continuous Renal Replacement Therapy Using a Novel Biomarker in Patients with Acute Kidney Injury after Cardiac Surgery-An Observational Trial. *International journal of molecular sciences*. 2024;25(20).
16. Benning L, Reineke M, Rodrigues CE, Kalble F, Speer C, Sommerer C, et al. Proenkephalin A 119-159 in Kidney Transplantation: A Novel Biomarker for Superior Tracking of Graft Function Trajectories. *Transpl Int*. 2025;38:14366.
17. Kieneker LM, Hartmann O, Struck J, Bergmann A, Gansevoort RT, Joosten MM, et al. Plasma Proenkephalin and Poor Long-Term Outcome in Renal Transplant Recipients. *Transplantation direct*. 2017;3(8):e190.
18. Ponce D, Balbi A. Acute kidney injury: risk factors and management challenges in developing countries. *Int J Nephrol Renovasc Dis*. 2016;9:193-200.
19. Hoste EA, Bagshaw SM, Bellomo R, Cely CM, Colman R, Cruz DN, et al. Epidemiology of acute kidney injury in critically ill patients: the multinational AKI-EPI study. *Intensive care medicine*. 2015;41(8):1411-23.
20. Kellum JA, Sileanu FE, Bihorac A, Hoste EA, Chawla LS. Recovery after Acute Kidney Injury. *American journal of respiratory and critical care medicine*. 2017;195(6):784-91.
21. Hall PS, Mitchell ED, Smith AF, Cairns DA, Messenger M, Hutchinson M, et al. The future for diagnostic tests of acute kidney injury in critical care: evidence synthesis, care pathway analysis and research prioritisation. *Health Technol Assess*. 2018;22(32):1-274.
22. Ostermann M, Zarbock A, Goldstein S, Kashani K, Macedo E, Murugan R, et al. Recommendations on Acute Kidney Injury Biomarkers From the Acute Disease Quality Initiative Consensus Conference: A Consensus Statement. *JAMA Netw Open*. 2020;3(10):e2019209.
23. Zarbock A, Forni L, Koyner JL, Gomez H, Pannu N, Ostermann M, et al. Preventing acute kidney injury and its longer-term impact in the critically ill. *Intensive care medicine*. 2025;51(7):1331-47.
24. Schaeffner E. Determining the Glomerular Filtration Rate-An Overview. *J Ren Nutr*. 2017;27(6):375-80.
25. Levey AS, Inker LA, Coresh J. GFR estimation: from physiology to public health. *American journal of kidney diseases : the official journal of the National Kidney Foundation*. 2014;63(5):820-34.
26. Shemesh O, Golbetz H, Kriss JP, Myers BD. Limitations of creatinine as a filtration marker in glomerulopathic patients. *Kidney international*.

1985;28(5):830-8.

27. Beunders R, van Groenendaal R, Leijte G, Kox M, Pickkers P. Proenkephalin Compared to Conventional Methods to Assess Kidney Function in Critically Ill Sepsis Patients. *Shock*. 2020.
28. Khorashadi M, Beunders R, Pickkers P, Legrand M. Proenkephalin: A New Biomarker for Glomerular Filtration Rate and Acute Kidney Injury. *Nephron*. 2020;1-7.
29. Patschan D, Erfurt S, Oess S, Lauxmann M, Patschan S, Ritter O, et al. Biomarker-Based Prediction of Survival and Recovery of Kidney Function in Acute Kidney Injury. *Kidney Blood Press Res*. 2023;48(1):124-34.
30. Dépret F, Hollinger A, Cariou A, Deye N, Vieillard-Baron A, Fournier MC, et al. Incidence and Outcome of Subclinical Acute Kidney Injury Using penKid in Critically Ill Patients. *American journal of respiratory and critical care medicine*. 2020;202(6):822-9.
31. Grossman A, Clement-Jones V. Opiate receptors: enkephalins and endorphins. *Clin Endocrinol Metab*. 1983;12(1):31-56.
32. Denning GM, Ackermann LW, Barna TJ, Armstrong JG, Stoll LL, Weintraub NL, et al. Proenkephalin expression and enkephalin release are widely observed in non-neuronal tissues. *Peptides*. 2008;29(1):83-92.
33. Sezen SF, Kenigs VA, Kapusta DR. Renal excretory responses produced by the delta opioid agonist, BW373U86, in conscious rats. *The Journal of pharmacology and experimental therapeutics*. 1998;287(1):238-45.
34. Grossman A, Besser GM, Milles JJ, Baylis PH. Inhibition of vasopressin release in man by an opiate peptide. *Lancet*. 1980;2(8204):1108-10.
35. Ernst A, Kohrle J, Bergmann A. Proenkephalin A 119-159, a stable proenkephalin A precursor fragment identified in human circulation. *Peptides*. 2006;27(7):1835-40.
36. penKid(R) SR. Instructions for Use (IFU), IFU-PEK-001 (EUen). March 2002.
37. Kim H, Hur M, Lee S, Marino R, Magrini L, Cardelli P, et al. Proenkephalin, Neutrophil Gelatinase-Associated Lipocalin, and Estimated Glomerular Filtration Rates in Patients With Sepsis. *Annals of laboratory medicine*. 2017;37(5):388-97.
38. Smeets NJL, Hartmann O, Schulte J, Schreuder MF, de Wildt SN. Proenkephalin A as a marker for glomerular filtration rate in critically ill children: validation against gold standard iohexol GFR measurements. *Clinical chemistry and laboratory medicine : CCLM / FESCC*. 2022;61(1):104-11.
39. KDIGO. KDIGO clinical practice guideline for acute kidney injury. *Kidney international*. 2012;2:1-138.
40. Zou C, Wang C, Lu L. Advances in the study of subclinical AKI biomarkers. *Front Physiol*. 2022;13:960059.
41. Dépret F, Polina A, Amzallag J, Fayolle-Pivot L, Coutrot M, Chaussard M, et al. PenKid measurement at admission is associated with outcome in severely ill burn patients. *Burns*. 2020.
42. Casalboni S, Valli G, Terlizzi F, Mastracchi M, Fidelio G, De Marco F, et al. 30 Days Mortality Prognostic Value of POCT Bio-Adrenomedullin and Proenkephalin in Patients with Sepsis in the Emergency Department. *Medicina (Kaunas)*. 2022;58(12).
43. Papasidero ID, Valli G, Marin D, Del Sasso A, De Magistris A, Cennamo E, et al. Utility of Measuring Circulating Bio-Adrenomedullin and Proenkephalin for 30-Day Mortality Risk Prediction in Patients with COVID-19 and Non-COVID-19 Interstitial Pneumonia in the Emergency Department. *Medicina (Kaunas)*. 2022;58(12).
44. Schiffli H, Lang SM. Current Approach to Successful Liberation from Renal Replacement Therapy in Critically Ill Patients with Severe Acute Kidney Injury: The Quest for Biomarkers Continues. *Mol Diagn Ther*. 2021;25(1):1-8.
45. Sigurjonsson J, Nyman U, Bjursten H, Grubb D. Renal elimination and eGFR prediction of proenkephalin. *Scandinavian journal of clinical and laboratory investigation*. 2025;85(4):281-6.
46. Lorenzin A, de Cal M, Perin N, Morisi N, Brendolan A, Lentini P, et al. Human Proenkephalin A 119-159 (penKid) in Extracorporeal Therapies: Ex vivo Sieving Coefficient, Diffusive Clearance, and Hemoadsorption Kinetics. *Blood Purif*. 2024;53(10):773-80.
47. Yarlagadda SG, Coca SG, Formica RN, Jr., Poggio ED, Parikh CR. Association between delayed graft function and allograft and patient survival: a systematic review and meta-analysis. *Nephrology, dialysis, transplantation : official publication of the European Dialysis and Transplant Association - European Renal Association*. 2009;24(3):1039-47.
48. Lin LC, Chuan MH, Liu JH, Liao HW, Ng LL, Magnusson M, et al. Proenkephalin as a biomarker correlates with acute kidney injury: a systematic review with meta-analysis and trial sequential analysis. *Critical care*. 2023;27(1):481.

Your notes

Disclaimer

This communication is for educational use only. It is not intended to guide or imply medical decision making.

sphingotest® penKid® is offered for research use only. “PenKid” represents the analyte Proenkephalin A 119-159.

For more information, please visit www.sphingotec.com

SphingoTec GmbH
Neuendorfstr. 15A
16761 Hennigsdorf, Germany
Phone: +49 3302 20565 0
Email: info@sphingotec.com
www.sphingotec.com